Don't Be Basic: Connecting laboratory testing to the field via evaluation of an advanced silicone foam dressing

Authors[†]: Laura Swoboda, DNP

Aligning foam dressing performance with clinical goals

Problem and Purpose: The development of effective and cost-efficient wound dressings is a significant area of research in wound care. Silicone bordered foam dressings have emerged as a popular choice due to their ability to maintain a moist wound environment, absorb exudate, and minimize trauma during dressing changes. This report on the translation of laboratory findings and clinical outcomes aims to describe the performance of a novel silicone bordered foam dressing in terms of its adhesive properties, fluid handling, and costeffectiveness.

Methods: A series of laboratory-based studies¹ were conducted to assess the performance characteristics of the novel dressing² compared to other common dressings types in both regular and sacral designs (n=10). Adhesive strength, friction, shear were measured using a ground steel sled to determine forces-dressing interactions. Fluid handling was evaluated by exposing the dressing in simulated wound exudate and measuring fluid handling capacity. The costeffectiveness of the dressing was assessed by comparing its price to other commercially available silicone bordered foam dressings.∞

Average Selling Price = (ASP)

Product Revenue Total number of units sold

∞ Based on Avg. Sales Price (ASP) from Clarivate – DRG data (decisionresourcesgroup.com)

Results: The novel silicone bordered foam dressing demonstrated excellent adhesive properties, adhering firmly to the skin simulants without over securement, which can be a concern for medical adhesive related skin injury and dermatitis. It exhibited a high absorbency capacity, effectively managing a wide range of exudate levels with appropriate moisture vapor transition rate to manage the dermal-dressing microclimate. Additionally, the dressing proved to be cost-effective, offering comparable performance to more expensive brands.

Name Point deflection k-----Strain ratio k-----

Test

Fluid Handling Capacity

Distance of nonconformity

Product Name ALLEVYN LIFE

OPTIFOAM GENTLE EX ConvaFoam MEPILEX BORDER FLEX HALYARD* ALLEVYN LIFE OPTIFOAM GENTLE EX ConvaFoam HALYARD*

Clinical Purpose/ Method Correlations Figure 5. Comparison of temperature at 175 minutes (Sacral Dressings) Error bars indicate a 95% confidence interval Characterizes The resulting cone-shaped Figure 1. Sacral dressings: Force after 30 seconds. Error bars indicate a 95% confidence interval. the effect of a deflection demonstrates the **Temperature at 175 minutes Average Force at Finish** ability of the dressing to perpendicular distribute that force to force applied to a Allevyn Life Sacrum wound dressing. surrounding areas, away 35.5 Optifoam Gentle EX 35.0 ntribute from the point of force. 34.5 ealing. Mepilex Border Sacrum 34.0 33.5 **9** 33.0 ConvaFoam Silicone 32.5 32.0 HALYARD* Bordered 31.5 Silicone Foam 31.0 Control nside Dressing (Indenter-Outside Dressing (Dressing-Support Surface Interface) Dressing Interface) Allevyn Lite ConvaFoam Mepilex Border HALYARD^{*} Border Silicone Foam Figure 6. Percentage of forces transmitted through the sacral dressings. Error bars indicate a 95% confidence interval. Coefficients of Friction: 6x6 Dressings Ratio based on The strain ration tells you i Figure 2. Sacral dressings: Average strain ratio. Error bars indicate a 95% confidence interval. the deflection eper a dressing can spread the Average Strain Ratio force of a boney prominence diameter and the 0.50 deflected materia like the sacrum and coccyx length. Closer to from a small area of focus at 0.850 1 indicates the the tip of the bone to wider 0.30 dressing is areas and less pressure 0.800 distributing the overall. This is another 0.20 0.750 applied force to reason why we use larger 0.700 0.10 dressings for PIP greater a wider area. than the area directly over 0.650 the boney prominence. Optifoam Gentle EX ConvaFoam Silicone Mepilex Border Flex HYH Bordered Allevyn Life 0.600 Silicone Foam 0.550 Static Friction Coefficient Kinetic Friction Coefficient Figure 7. Comparison of coefficients of friction for sacral dressings 0.500 Error bars indicate a 95% confidence interval (n=15) Mepilex Border HALYARD* Bordered Control Allevvn Life ConvaFoam FX Silicone Sacrum Silicone Foam Coefficients of Friction: Sacral Dressings Dual-use (PIP and wound Measuring and Figure 3. Comparison of fluid handling capacity of small dressings (4x4 inches). und 0.50 comparing treatment) dressings must Error bars indicate the standard deviation. the Moisture handle wound drainage and 0.40 Vapor Loss, perspiration while preventing **Fluid Handling Capacity** Fluid Absorbed, significant contamination 0.30 2.5 and total fluid from urine, feces, and 0.20 environmental sources. handling capacity. HALYARD* Bordered 2.0 Fluid can be passed out Silicone (4x4) 0.10 through the dressing, or ConvaFoam Silicone the dressing can retain the 0.00 <u>∼</u> 1.5 Allevyn Life fluid within. Silicone Foam Mepilex Border Flex (4x4)Static Friction Coefficient
Kinetic Friction Coefficient Optifoam Gentle EX adds (4×4 Figure 8. Comparison of pressure map images Allevyn Life Sacrum Opti foam Gentle EX Allevyn Life (4x4) Control clinical 5 10 15 20 25 30 35 40 45 50 5 10 15 20 25 30 35 40 45 50 5 10 15 20 25 30 35 40 45 50 sing's Fluid Absorbed Moisture Vapor Loss **Total Fluid Handling** ressure Capacity sing and n of Figure 4. Dressing non-conformity distance comparison. Error bars indicate a 95% confidence interval Dressing Distance of nonconformity non-conformity is the distance between the Distance of Non-Conformity intergluteal cleft and the is measured dressing. Gaps in dressings at the intergluteal cleft of a heated adherence leave entrances 7.0 for contaminants and leakage mannequin before and after for wound exudate. Gaps can 6.0 HALYARD^{*} Bordered also prevent the dressing application of a ConvsFoam Silicone Mepilex Border Sacrun Silicone Foam 5.0 5 10 15 20 25 30 35 40 45 50 5 10 15 20 25 30 35 40 45 50 10 15 20 25 30 35 40 45 5 from adequately managing clinically relevant pressure forces. load. 3 4.0 3.0 2.0 1.0 Silicone Foam, Sacrum Silicone, Sacral 5.0 - 50 O Minutes 15 Minutes 30 Minutes Average SKU **Conclusion:** The novel silicone bordered foam dressing² shows promising results in terms of its adhesive properties, absorbency 66801067 capacity, and cost-effectiveness. These findings suggest that this dressing may be a valuable addition to the wound care arsenal, MSCEX44EP 423253 providing a reliable and affordable option for wound management. Further clinical studies are warranted to evaluate the dressing's 595300 performance in real-world settings and assess its impact on wound healing outcomes. 50085 66801306 [†] EC Service, Inc. – Evan Call, Kasey Call, Sandra Guzman, Marianne Russon MSCEX77EP 1. EC Service, Inc., Centerville, UT, USA 2. O&M Halyard, Alpharetta, GA, USA

- MEPILEX BORDER SACRUM
- 423256 282055 50087





Test Name	Purpose/ Method	Clinical Correlations
Temperature	Temperature was measured at the dressing-mannequin interface (inside dressing) and dressing-support surface interface (outside dressing) at 175 minutes after dressing was applied to a heated mannequin with a clinically relevant load.	Human skin maintai a temperature just le than core body temperature. Elevati and depressions in temperature can cor to delayed wound he tissue breakdown, a increase risk of wou occurrence.
Shear	Simulated boney prominence via a ground steel sled with an angled face slides over the dressing's surface.	Shear can contribute pressure injury form as it impacts the dee tissue levels. It can a impact perfusion and wound healing.
<section-header></section-header>	Static and Kinetic co-efficients of friction are measured alongside shear.	Friction contributes to pressure injury forms especially in the more superficial tissues. It can also disturb would healing.
Pressure mapping	XSENSOR pressure mapping technology is used to evaluate sacral mannequin w/wo product.	Pressure mapping a another layer of measurement to the evaluation of a dress ability to disperse pr throughout the dress prevent the formatio pressure injuries.
Cost analysis	Publicly available pricing via a USA disruptor was utilized to compare pricing of dressings.	It may be incorrectly assumed that dressi performance in correct to the price of that dressing. Comparing pricing amongst dressings alongside their performance characteristics provis of objective guidance clinicians when task with making inventor decisions.









