

The Clinical Challenge in Diabetic Foot Ulcers



MAP Wound Matrix, a Revolutionizing Scaffold

Microporous Annealed Particle Post-Mohs wound (MAP) Wound Matrix is a 'smart' flowable and resorbable wound matrix with a unique porous microenvironment that enables rapid tissue integration without eliciting a foreign body response. The MAP Wound Matrix enables, with a single application, a new treatment for wounds of any shape without the use of cells, biologics or growth factors.

MAP Wound Matrix is delivered to the wound as a flowable (1) Flowable 100 μ m concentrated hydrogel microsphere solution, and then,



building blocks

upon white light exposure, is transitioned in situ to a robust scaffold. The resultant scaffold displays a functional open pore structure with microscale porosity (pore size: 10 - 50 µm) to allow free movement of cells to build new tissue prior to enzymatically mediated material resorption.



Objective of the Study

The purpose of this study was to evaluate the efficacy of the MAP Wound Matrix, a novel biomaterial, as a single application treatment to promote diabetic wound healing using a diabetic porcine model.

35-Day Full-thickness Wound Healing Study of MAP Wound Matrix in Diabetic Yucatan Miniature Swine

Donald R. Griffin PhD^{1,2}, Stephanie Deshayes PhD¹, and Westbrook M. Weaver PhD¹.

¹Tempo Therapeutics, Inc., San Diego, CA, ²Biomedical Engineering, Chemical Engineering, University of Virginia, Charlottesville, VA.

Introduction

An estimated 20% of DFU patients require lower limb amputation.

Patients with DFU have a 5-year patient mortality rate of 74%.

These statistics outline a clear need for innovative wound care materials that can accelerate healing and improve outcomes.



References

2. D.R. Griffin, et al, Nature Materials, 2021 20: 560–569.

Methods

Twelve full-thickness wounds (2 x 2 cm) were created on the back of each alloxan-induced diabetic Yucatan miniature swine. Wounds were treated with either MAP Wound Matrix (referred as MAP), Aquaphor (standard of care), or Oasis SIS (comparator). Animals were followed for up to 35 days. This study measured the histological wound healing characteristics as primary endpoints including re-epithelialization, inflammation, vascularization, and device degradation. The study also measured clinical wound observations as secondary endpoints, such as skin reactions, and wound closure.



The MAP Wound Matrix holds promise as an alternative therapeutic option for DFUs, including the promotion of multiple superior wound healing outcomes relative to current clinical options. Importantly, the scaffold's ability to act as a single treatment option that reduces the need for repeated clinician intervention suggests that it may improve patient outcomes and reduce healthcare burden. Further clinical studies will validate its efficacy and safety in patients.

