

# The role of nitric oxide in the killing & prevention of surface-associated bacterial communities by a prototype nitric oxide-generating wound dressing

Daniel Metcalf PhD<sup>1\*</sup>, Alan Horner PhD<sup>1</sup>, Nicholas Boote PhD<sup>1</sup>

<sup>1</sup> Convatec Ltd, Oxford, UK

## Introduction

- Hard-to-heal wounds, such as diabetic foot ulcers (DFUs), are often compromised by microorganisms that contribute to chronicity and infection risk, particularly in high-risk patients such as those with diabetes<sup>1</sup>
- Nitric oxide (NO) is a potent antimicrobial<sup>2</sup> and antibiofilm<sup>3</sup> gas produced by the mammalian innate immune response to microorganisms, yet it has unrealised potential in wound care<sup>4</sup>
- A novel NO-generating dressing (NOGD) showed superiority over standard of care (SoC) in a randomised controlled trial (RCT) in DFUs<sup>5</sup>, and has demonstrated antibiofilm activity *in vitro*<sup>6</sup>
- NO is generated in NOGD when two component parts are placed together: a Carrier Layer (CL) containing aqueous nitrite and a low pH Absorbent Protective Layer (APL; which acidifies nitrite to NO)

Objective

The aim of these *in vitro* studies was to examine the antibiofilm activity of the NO-generating technology within NOGD

## Methods

- Wound-associated pathogens were used in biofilm kill and prevention assays, while varying nitrite concentration and treatment times
- Initial prototype test dressings comprised APLs plus CLs containing varying concentrations of nitrite
- **Nitrite (NO precursor) concentration:**
  - Biofilm kill: biofilms of methicillin-resistant *Staphylococcus aureus* (NCTC 12493; MRSA) were grown on 25 mm dia. nitrocellulose filters on agar for 24 hours, as described elsewhere<sup>6</sup>, before prototype and controls (no dressing, and dressings without nitrite) were applied
  - Biofilm prevention: 1x10<sup>5</sup> CFU/mL planktonic MRSA-inoculated filters were immediately treated with prototypes and controls for 24 hours
    - In both assays, surviving bacterial cells were enumerated in quadruplicate on agar
- **Treatment duration:**
  - The final NOGD design containing a 1 M nitrite CL was used to regularly measure biofilm survival or formation over a 24-hour period, as described above and elsewhere<sup>6</sup>

## Results

### Nitrite concentration:

- MRSA biofilm was reduced by 3 log<sub>10</sub> in 24 hours by prototypes with CLs containing 0.2 M nitrite, and was eradicated by prototypes with CLs containing 0.5 M nitrite (Fig 1)
- MRSA biofilm formation was not prevented by prototypes with CLs containing 0.1 M nitrite, but was completely prevented by prototypes with CLs containing 0.2 M nitrite, after 24 hours (Fig 2)
  - This data suggests the NO generated by 0.5 M and 1 M nitrite in the CL is sufficient to kill and prevent MRSA biofilm *in vitro*




Figure 1. MRSA biofilm kill by dressing prototypes containing nitrite at varying concentrations (N=4). Placebo=CL with no nitrite




Figure 2. MRSA biofilm prevention by dressing prototypes containing nitrite at varying concentrations (n=4). Placebo=CL with no nitrite

### NOGD treatment duration:

- Using NOGD with CL containing 1 M nitrite, MRSA biofilm was reduced by 3 log<sub>10</sub> after 2 hours, by >7 log<sub>10</sub> after 4 hours, and completely eradicated after 6 hours (Fig 3)
- Biofilm formation was completely prevented after 6 hours (Fig 4)
  - This data confirms the NOGD can completely kill MRSA biofilm in 6 hours, and prevent its formation in 6 hours, in an NO-dose dependent manner, *in vitro*




Figure 3. MRSA biofilm kill by NOGD after various treatment durations. (■) control; (■) placebo with no nitrite; (■) NOGD (N=4)




Figure 4. MRSA biofilm prevention by NOGD after various treatment durations. (■) control; (■) placebo with no nitrite; (■) NOGD (N=4)

## Discussion

- Biofilm kill and prevention models show that varying the concentration of nitrite used to generate NO, and varying the duration of treatment with the final NOGD design, resulted in a NO dose-response effect on MRSA biofilm kill and prevention
- Nitric oxide, a natural antimicrobial molecule of our innate immune system with intrinsic antibiofilm activity, presents a low risk of antimicrobial resistance<sup>7</sup> and biofilm tolerance<sup>4</sup>, unlike many other antimicrobial agents employed in wound care today
- These studies build on significantly superior DFU clinical outcomes<sup>5</sup> and early *in vitro* performance<sup>6</sup> of this novel dressing technology

Conclusion

Nitric oxide within an advanced wound dressing has high potential as an antibiofilm agent in the treatment of hard-to-heal wounds

Acknowledgements: The laboratory work presented here was conducted by Richard Waite and Robert P Allaker, Queen Mary University of London

1. Armstrong et al. Diabetic foot ulcers and their recurrence. N Engl J Med 2017; 376: 2367-237  
2. Bogdan C. Nitric oxide and the immune response. Nature Immunol 2001; 2: 907-916  
3. Barraud et al. Nitric oxide-mediated dispersal in single- and multi-species biofilms of clinically and industrially relevant microorganisms. Microb Biotechnol 2009; 2: 370-378  
4. Roberts et al. Harnessing the power of our immune system: the antimicrobial and antibiofilm properties of nitric oxide. Microorganisms 2024; 12: 2543  
5. Edmonds et al. Multicenter, randomized controlled, observer-blinded study of a nitric oxide generating treatment in foot ulcers of patients with diabetes-ProNOx1 study. Wound Repair Regen 2018; 26: 228-237  
6. Waite et al. Activity of a nitric oxide-generating wound treatment system against wound pathogen biofilms. Int J Antimicrob Agents 2018; 52: 338-343  
7. Privett et al. Examination of bacterial resistance to exogenous nitric oxide. Nitric Oxide 2012; 26: 169-173