

كليسة العلسوم والتقنسيات فساس

+.YEU.I+ I +C.OO.IEI A +OIEREFEI - X.O FACULTÉ DES SCIENCES ET TECHNIQUES DE FÈS

3D-Printed Antimicrobial Wound Dressing with Argan Nutshell Lignin & Essential Oils for Advanced Healing

Khaoula Sebbar (1), Masoud Adhami (2), Eneko Larraneta (2)*, Nada Al Moudani (1), Lotfi Aarab(1), Saad Koraichi Ibnsouda (1), Soumya Elabed (1)**

(1) Laboratory of Microbial Biotechnology and Bioactive Molecules, Faculty of Sciences and Technologies, Sidi Mohamed Ben Abdellah University of Fez BP 2202, Morocco.

(2) Belfast School of Pharmacy, Queen's University, 97 Lisburn Road, Belfast BT9 7BL, United Kingdom





→ Materials :

Solubilize

MATERIALS & METHODS

Antimicrobial Agents: MPEO and ROEO,

- Biopolymer: Lignin Extracted from Agan,

- Synthetic polymer: Eudragit RS PO.


CHRONIC WOUNDS







- Natural resources and Agro-waste in developing countries like Morocco are poorly



2. Design of Cost-Effective, Tailor-Made Meshes via 3D Printing

3. Replacement of Synthetic Antibiotics with Natural Agents:

4 Integration of Fudragit as a Smart Polymeric Matrix: imploy Eudragit (EUD) for its cont



→ 3D Printing Formulations:

Table 1: Composition of the formulations that go into making each group of the 3D printed wound

ell Culture :		Composition (N)						
B22°C (C)	Fermulations	Polymers				Solvent	Essential O	
		ANS-UG	EUD ES PO	res	Soluplus	Acetone	MPEO	*
	Ave EUO		ASN	5%	10%	35N		
	EU0-469E0		75N	SN	10%	35N	30%	
	EU0-9000		75N	5%	10%	JSN		
	TWARRO	40%	JSN	SN	20%	JSN		
17°C 55°CO,	THEALTONIMACO	ADN	JSN	SN	10%	JIN	20%	
Fibroblast culture	110/200-4000	40%	JSN	5%	10%	35N		

→ Step by Step Process : → Characterizations of the 3D printed Wound Dressings:




1. Morphological Characterizations 2. Chemical Characterization 3. Mechanical Characterization 4. Thermal Characterizations 5. Contact Anale Measuremen Wound Dressings

> 6. Antioxidant Activity 7. Swelling and Degradation Behaviors 8. Antibacterial Properties


> > 10. In Vitro Cytotoxicity Test

RESULTS & DISCUSSIONS



pure EUD, MPEO, and ROEO ensuing 3DP samples of EUD control, LIG/EUD, EUD+EOs, and

Figure 5: TGA thermograms of pure LIG. pure EUD. MPEO, and FUD control LIG/FUD FUD+FOs and LIG/FUD+FOs




Printed Wound Dressings Infused with MPEO or ROEO, Starting from a 50 mg/L DPPH Concentration.

Figure 3: ATR-FTIR spectra of EUD powder with or w/o LIG, 3D-printed wound dressings containing MPEO or ROEO.



Figure 6: Assessment of Eudragit RS PO and Argan Shell Lignin Composite Wound Dressings, with and without MPEO or ROFO: (a) Contact Angle Visualization and (b) Quantitative



(a), E. coli (b), S. aureus (c), and P. aeruginosa (d) through Optical Density (OD) Assay

CONCLUSIONS

Table 2: The tensile and flexural mechanical characteristics

Figure 7: Swelling Kinetics and Degradation Behavior of the

3D-Printed Wound Dressings.

nvolving the different tested wound dressings.

- First-time use of Moroccan Argan Nutshell Lignin (ANS-LIG) was Successful. by valorizing local agro-waste as a 3D-printed biopolymer for biomedical use.
- Strong Antibacterial Performance using Essential Oils, especially LIG/EUD+MPEO Wound Dressing, which drastically reduced S. gureus (6%) and E. coli (7.5%) viability. proving its potency against wound pathogens.
- Excellent Biocompatibility and Dermal Cell Proliferation, since all of the 3Dprinted Patches enhanced fibroblast viability, especially LIG/EUD+ROEO with up to 225%, confirming safe and potential regenerative interaction with skin cells.

AKNOWLEGEMENTS

CONTACT

Khaoula SEBBAR **PhD Student**

FST-Fez. MOROCCO

of Control and 3D-Printed EUD-Based Scaffolds Enriched with